




LinkSymphoKnee
Built on TRUST



“LINK is no longer just a revision knee company for me. The new LINK knee is my go-to for primary cruciate retaining total knee arthroplasty. The narrow nature of the implant and design of the trochlea groove gives me an implant that fits all my patients.”


– Brett Levine MD Midwest Orthopaedics at RUSH¹

¹Dr. Brett Levine performed the first global implantation of the LSK cruciate retaining design in Chicago, Illinois in August of 2021.



LinkSymphoKnee

Built on TRUST



Dual Radius of Curvature Designed to Optimize Kinematics

High Flex Recess for Patellar Ligament


Flatter Mediolateral Edge at Posterior Facet for Up Down Femoral Compatibility

Forgiving Anterior Profile

3° Posterior Slope

Posterior and Central Dovetail Locking Feature

E-Dur Vitamin E Infused Polyethylene¹



¹Crosslinking w/ a 75 kGy gamma radiation dose in air then EtO sterilized.